



The „Quickstep“ to shoe lasts



Push in. Press in.
Press on... finished!



From Germany. And developed from practice.

Consulting and bureaucratic work are increasingly becoming the focus of everyday work for orthopaedic shoe technicians. There is less and less time for the craftwork that really matters. Nevertheless, the level of quality must remain high. This is one of the reasons the Kuickbox was developed: to save time while producing excellent impressions.

The Kuickbox was developed from orthopaedic shoe technicians to make easier the central step for precise last producing: the foot impression. Previous procedures are partially too complicated or time consuming – and to that related with high acquisition costs. With the Kuickbox it works straightforward faster, easier and in TOP-Quality.



The Quickbox: Foot impression quickly made!

Quickbox is the new device for fast production of lasts for custom made shoes or orthopaedic lasts. Detailed, clean and cost-effective. For low shoes, ankle or leg lasts.

The advantages at a glance

- super fast and accurate moulding
- detailed depiction of all foot and leg details
- correction of abnormal foot positioning during moulding
- very easy to use
- no drying times
- minor reworking

No subsequent casting on of the toe-caps

- saves a lot of time

Satisfied patients

- clean work on the foot
- innovative process
- impressive design

Manageable costs

Easy to transport

Integrated impression control function

Delivery program:

- 770 043 01** Quickbox Sz. 1 low shoe last
- 770 043 02** Quickbox Sz. 2 foot last
- 770 043 11** Foam for refilling Sz. 1 low shoe
- 770 043 12** Foam for refilling Sz. 2 foot
- 770 043 15** Functional wedges and height adjustment plates

Made in Germany



Götz Service GmbH · Zillenhardtstr. 2/2 · DE-73037 Göppingen-Voralb
phone +49(0)7161.6102-0 · fax +49(0)7161.6102-990
info@goetz-service.com · www.goetz-service.com